

Selecting by marquee1

Adjust the preferences for object and border colors, units of measure, and so on.
Marquee selection allows you to more easily select several objects at a time.

1. Click the **Plate Builder Tool** button to start a Plate Builder session.
2. On the Plate Builder toolbar, click **Select Object**.
3. To select several objects at one time, hold the mouse button and drag a selection rectangle (a marquee) completely around the objects that you want to select.
Note: Plate Builder selects only objects that are fully inside the marquee.
A selected object appears highlighted with the color that you specified.
4. To add an object, hold down Shift, and drag a marquee around the objects.
5. To remove objects from a selection group, hold down Shift, and drag a marquee completely around the objects that you want to remove from the group.