

_SvIncludePlusContent-Calibrating the color monitors

You must calibrate your monitors when the calibration expiry interval that you defined expires to ensure accurate color viewing. Monitors must be turned on for at least 10 minutes before they can be calibrated.

Note: The MV calibration tool checks for stability. If a sufficient warm-up time has not been achieved, calibration will not be successful.

1. At the Dock, click the **Matchprint Virtual folder** icon.
2. Select **MVCalibrator**.
During calibration, Matchprint Virtual checks the monitor's color and luminance and compares it to a reference profile for that monitor type. If the monitor cannot achieve the specified luminance or achieve excellent color match, the calibration fails.